
Levenhuk MED PRO 600 Fluo Microscope

- EN User Manual
- BG Ръководство за потребителя
- CZ Návod k použití
- DE Bedienungsanleitung
- ES Guía del usuario
- HU Használati útmutató
- IT Guida all'utilizzo
- PL Instrukcja obsługi
- PT Manual do usuário
- RU Инструкция по эксплуатации
- TR Kullanım kılavuzu

levenhuk[®]
Zoom&Joy

1

EN

1	Trinocular head	15	Knob to adjust halogen lamp position in vertical direction
2	Luminous flux switch (splitter)	16	Collector adjustment knob
3	Trinocular head lock knob	17	Lamp filament adjustment knob
4	Field diaphragm alignment knob	18	Power on/off button
5	Field diaphragm opening adjustment knob	19	Fine focusing knob
6	Aperture diaphragm opening adjustment knob	20	Knob to move slide in cross direction
7	Collector adjustment knob	21	Knob to move slide lengthwise
8	Mercury-filled lamp flashlight	22	Condenser
9	Cap clamping screw	23	Field diaphragm opening adjustment ring
10	Knob to adjust mercury-filled lamp position in vertical direction	24	Condenser clamping screw
11	Knob to adjust mercury-filled lamp position in horizontal direction	25	UV protection screen
12	Microscope stand	26	UV protection screen clamping screw
13	Halogen lamp flashlight	27	Switching ring
14	Knob to adjust halogen lamp position in horizontal direction	28	Fluorescent illuminator clamping screw

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Microscope description and operation

Application

The microscope is designed for diagnostic testing, including by immunofluorescence method, in clinical, microbiological, pathoanatomical, and other laboratories in medical institutions. Besides, it may also be used in veterinary science, crop growing, bioengineering, pharmaceutical industry, for expertise in the sphere of criminalistics, state epidemiological surveillance, environmental protection. The microscope is used to study stained and unstained slides in the form of smears and microsections in transmitted light.

In the luminescence light, the microscope makes it possible to detect dangerous bacterial and virus infections when observing objects stained by Auramine, Acridine orange, FITC etc.

Being properly operated, the microscope is safe for health, life, property of the consumer and for the environment. The microscope stand is of anti-vibration design. The microscope is made to operate at ambient air temperature from +10 to +35°C (50 to 95°F) and relative humidity of 80% max. The oil-immersion objective lens shall be operated indoors at ambient air temperature from +15 to +25°C (59 to 77°F).

Microscope design and operation principle

TO PREVENT MICROSCOPE BREAKDOWN, BEFORE STUDIES, CAREFULLY REVIEW THE RULES OF HANDLING AND PROCEDURE OF WORKING WITH THE MICROSCOPE SPECIFIED IN THIS OPERATION MANUAL.

Fluorescent microscope operating principle is based on using a phenomenon of fluorescence (luminescence) of observed objects caused by rays of light having a specific spectrum. To excite fluorescence, objects are illuminated from the top through the objective lens, with a mercury-filled lamp used as a source of such light. The luminous flux required to excite fluorescence is separated from the total radiation of the mercury-filled lamp with filters, conventionally referred to as the excitation filters.

To guide the luminous flux to the objective, a beam splitter with a special interference coating is used that mostly reflects excitation light and transmits object fluorescence light. The excitation filter, the beam splitter and the cut-off filter (used to absorb the residual excitation radiation) are combined in a single beam splitting unit. A kit of five beam splitting units is mounted on a turret having a free socket for operating in transmitted light.

The optic system enabling object study in fluorescence light is made in the form of a detachable illuminator installed on the microscope stand. The microscope stand enables observation of objects illuminated with transmitted light.

Description and operation of components

Microscope stand

The microscope stand (fig. 1, 12) has ergonomic and stable shape, and is made of metal.

There is a double-stage focusing mechanism for vertical movement of the bracket (fig. 2, 3) with a coordinate stage (fig. 2, 10) and a revolving nosepiece for objective lens fixation (fig. 2, 7) on the stand. A fluorescent illuminator (fig. 2, 15) is installed on top of the stand and is fixed with a screw (fig. 1, 28). The microscope stand base contains the transmitted light illumination system and halogen lamp power source of 12V/30W. There is a socket for connection of a power cable on the rear surface of the base on the left.

The power source is built into the microscope stand base. Power on/off button (fig. 1, 18) energizes the halogen lamp installed in the flashlight (fig. 1, 13). Power is off in "O" position. Halogen lamp filament is adjusted by a knob (fig. 1, 17).

There is an iris field diaphragm, whose opening is adjusted by a ring (fig. 1, 23), located on the upper surface of the microscope base under the condenser (fig. 1, 22). A decorative base (fig. 2, 11) is put on the microscope stand base.

Revolving nosepiece

A six-position revolving nosepiece provides for objective lens (fig. 2, 7) setting into working parfocal position. The revolving nosepiece is inclined towards the microscope stand to provide space for installation and replacement of the examined slides.

Objective lenses are replaced by rotation of a corrugated ring (fig. 1, 27) of the revolving nosepiece to the fixed position.

Focusing mechanism

The focusing mechanism is designed for vertical displacement of the stage (fig. 2, 10) when the microscope is focused for sharp object image. The range of stage movement along height is 25mm. Vertical displacement of the stage is performed by coaxial knobs (fig. 2, 12 and 13) located on the left side of the microscope stand. The fine focusing mechanism knob (fig. 2, 12) has a scale with division value of 2µm. Behind the knob (fig. 2, 13) there is a ring (fig. 2, 14), designed to adjust ease of movement during coarse focusing. There is a fine focusing mechanism knob (fig. 1, 19) on the right side of the stand (fig. 1, 12).

Stage

A stage (fig. 2, 10) is equipped with a mechanism of coordinate object displacement in the horizontal plane in two mutually perpendicular directions. Stage and slide holder design (fig. 2, 8) provides for the possibility to install two slides and move them by 85mm in cross direction and 50mm lengthwise. Displacement is controlled by low-set coaxial knobs from the right side of the stand. With the help of the knob, the object is moved in the cross direction (fig. 1, 20) and lengthwise (fig. 1, 21). The division value is 1mm, the nonius division value is 0.1mm. The object is fixed on the stage surface between the holder and the clamp (fig. 2, 9) of the slide holder (fig. 2, 8). To install the object, the clamp (fig. 2, 9) is taken aside. The stage surface has a solid coating resistant to disinfection and wear. Stage dimensions are 180x160mm.

Transmitted light illumination system

The illumination system of the microscope is critical to obtain a contrast and evenly illuminated image of the objects under the microscope. The illumination system built into the microscope stand base (fig. 1, 12) is arranged according to the Köhler's principle in its classic version. The halogen lamp flashlight is installed on the rear wall of the base and is fixed with a screw (fig. 4, 2). The iris field diaphragm node is located on the stand base under the condenser (fig. 1, 22), the diaphragm opening is adjusted with a ring (fig. 1, 23).

The condenser (fig. 1, 22) is used to focus the field diaphragm image in the slide plane.

The illuminator is turned on with the help of a switch (fig. 1, 18) in position "I". The lamp brightness may be varied by rotation of the lamp filament adjustment knob (fig. 1, 17). The lamp is energized through a power supply source built into the microscope stand base.


The halogen lamp installation in the flashlight is shown in fig. 4. To access the lamp holder, it is necessary to unscrew the screw (fig. 1, 28) and to extend the cover (fig. 4, 5). To install the cover (fig. 4, 5) on the flashlight, it is necessary to take fixators (fig. 4, 6) behind the flashlight casing.

Bright-field condenser

An Abbe condenser (fig. 1, 22) is included in the microscope package for operation in the bright field. The condenser is installed into a bracket (fig. 3, 3) under the microscope stage and is fixed with a screw (fig. 1, 24).


An iris aperture diaphragm, diameter of which is adjusted with a knob (fig. 3, 1) changes aperture of the cone of rays illuminating the slide. A scale is applied onto the condenser frame, which makes it possible to reproduce illumination conditions selected for every objective lens – knob positions (fig. 3, 1).

There is an option to exclude the condenser frontal objective lens from the path of rays with the knob (fig. 3, 4) when working with low magnification objective lenses. Screws (fig. 3, 2) serve to align the field diaphragm image by displacement of the condenser in the horizontal plane. Condenser displacement along the optical axis of the microscope when focusing the field diaphragm image is carried out with the knob (fig. 3, 5).

- 1 Condenser aperture diaphragm opening adjustment knob
- 2 Condenser alignment screws
- 3 Condenser bracket
- 4 Condenser frontal objective lens switching knob
- 5 Condenser vertical displacement knob

Fig. 3: Condenser

- 1 Cap clamping seat
- 2 Flashlight clamping knob
- 3 Halogen lamp holder
- 4 Halogen lamp
- 5 Flashlight detachable cap
- 6 Fixators
- 7 Cap clamping screw

Fig. 4: Halogen lamp flashlight

Impinging light illumination system

The impinging light illumination system is made in the form of a detachable module – a fluorescent illuminator (fig. 2, 15). The module is installed with a lower flange into the microscope stand seat and is fixed with a screw (fig. 4, 7). A mercury-filled lamp flashlight is fixed on the illuminator (fig. 1, 8).

The fluorescent illuminator (fig. 2, 15) comprises a six-seat turret with 5 spectral beam splitting units and a free seat for transmitted light. Seats are numbered and provided with captions that contain information about spectral characteristics of filters and a dichroic mirror (a beam splitter) given in table.

	Excitation filter	Dichroic mirror	Cut-off filter	Unit designation
No.1	Free seat			
No.2	510–548	570	585–700	G
No.3	455–495	500	505–555	B
No.4	410–440	455	475	BV
No.5	380–440	435	450	V
No.6	330–370	405	425	U

Marking of beam splitting units complies with color of beam of rays exciting fluorescence of examined objects.

For example, when a carriage is in "G" position, a green spectral region of 510–560nm is identified from the total radiation flux of the mercury-filled lamp, and when it is in "B" position – a spectral region of 450–490nm (blue) is identified.

The illuminator comprises field and aperture diaphragms (FD and AD). The field diaphragm is equipped with an alignment device (fig. 1, 4), the position adjustment knobs are located on the illuminator casing on the right and left. The knob (fig. 1, 5) adjusts field diaphragm opening, the knob (fig. 1, 6) adjusts aperture diaphragm opening. To change the diaphragm dimensions, knobs (fig. 1, 5 and 6) are pushed into the casing. A plug to intercept light that is controlled by a knob is installed closer to the flashlight in the illuminator casing.

Mercury-filled lamp flashlight

The mercury-filled lamp flashlight (fig. 1, 8) is fixed on the illuminator with a bayonet ring as the adjusting ring (fig. 5, 7) and the fixator (fig. 5, 8) are brought closer to the illuminator end.

WARNING! PRIOR TO REMOVAL OF THE FLASHLIGHT FROM THE HEAD CASING, IT IS NECESSARY TO DISCONNECT THE MERCURY-FILLED LAMP POWER SUPPLY UNIT FROM THE GRID!

The mercury-filled lamp is aligned by knobs (fig. 1, 10 and 11). Knob 10 serves to move the holder with the lamp in vertical direction, and knob 11 – to move it in horizontal direction. The flashlight detachable cap (fig. 5, 1) is fixed with a screw (fig. 1, 9), on the inner side of which there is a mercury-filled lamp holder. The mercury-filled lamp (fig. 5, 4) is installed into bushings (fig. 5, 2 and 5) and is fixed by screws (fig. 5, 3 and 6).

WARNING! FOR MICROSCOPE TRANSPORTATION, REMOVE THE MERCURY-FILLED LAMP FROM THE FLASHLIGHT.

Inside the flashlight, there is a collector projecting the mercury-filled lamp discharge arc image into the exit pupil of the objective installed in the path of rays. Knob 7 (fig. 1) adjusts collector position along the illuminator axis.



Fig. 5: Mercury-filled lamp flashlight

- 1 Amperemeter
- 2 Lamp ignition button
- 3 Mercury-filled lamp operation-time counter
- 4 Power on/off button

Fig. 6: Mercury-filled lamp power supply unit

Trinocular head

The trinocular head (fig. 1, 1) is designed to work with objective lenses, the optical length of tubes in which is "infinity". The head enables binocular observation with eyepieces (fig. 2, 5) and image output through a vertical tube (fig. 2, 3) for fixation. The head enables positioning of eyepiece tubes in accordance with the required eyepiece spacing (observer's eye base). The distance between axes of eyepieces (fig. 2, 5) is adjusted by rotation of eyepiece tubes in the range from 50 to 75mm. The ring (fig. 2, 4) in the left eyepiece tube performs dioptric adjustment of eyepiece position (fig. 2, 5) for ± 5 diopters.

Luminous flux is directed into a vertical tube by switching the knob (fig. 1, 2). The knob can be switched into three positions, providing for three options of beam splitting: observation only, observation and documentation and only documentation. The head is installed into the seat of the fluorescent illuminator and is fixed by the fixator (fig. 1, 3).

Objective lenses

All objective lenses (fig. 2, 7) included into the scope of supply are designed for infinite optical length of the tube and have semi-apochromatic correction. Parfocal height of objective lenses is 45mm.

Linear magnification and numerical aperture are engraved on the casing of each objective lens. There is also a color marking compliant with magnification and information on a cover slip. The package includes objective lenses for working with the slides protected by cover slips and also objective lenses not requiring a slide to be protected by a cover slip.

Specifications of objective lenses

Correction type	Linear magnification and numerical aperture	System	Linear field in space of objects with eyepiece 10x/22, μm	Total microscope magnification with 10x eyepiece, x
Plan fluorite (semi-apochromatic)	4x/0.15	Dry	550	40
Plan fluorite (semi-apochromatic)	10x/0.35	Dry	220	100
Plan fluorite (semi-apochromatic)	20x/0.60	Dry	110	200
Plan fluorite (semi-apochromatic)	40x/0.75	Dry	55	400
Plan fluorite (semi-apochromatic)	100x/0.90 *	Dry	22	1000
Plan fluorite (semi-apochromatic)	100x/1.25 *	Oil	22	1000

* Not included into standard package

" ∞ /-" inscription on the objective lens means that the objective lens may handle slides both with a cover slip or without it. 40x and 100x magnification objective lenses are equipped with elastic frames preventing objects and frontal objective lenses against damage when focusing on object surface.

WARNING! IF OBJECTIVE LENSES ARE DAMAGED, THEY SHOULD BE REPAIRED IN THE MANUFACTURER'S SERVICE CENTER.

Eyepieces

The microscope package comprises two wide-field eyepieces with 10x magnification and 22mm linear field in the image plane.

Microscope use

Operation limits

The microscope should be used in premises where pushes and vibrations are hardly sensed, with no sources of intensive external exposure, i.e. sources of electromagnetic radiation. There shall be no excessive dust, acid, alkali vapors and other chemically active substances in the premises. The microscope should not be operated in brightly lit premises.

The microscope is designed for operation under the conditions of moderate and cold climate in laboratory premises at air

temperature from +10 to +35°C (50 to 95°F) and upper value of relative air humidity of 80% max.

Microscope unpacking

Unpack the microscope carefully and install it on a smooth surface. Check the microscope package contents. Visually inspect all elements included in the scope of supply, identify their purpose, make sure there are no damages and start assembling.

Microscope preparation for operation

Installation of modular units

- Install the decorative base on the microscope stand base (supply with the installed decorative base is possible).
- Install the halogen lamp flashlight (fig. 1, 13) on the microscope stand base and fix with a clamping knob (fig. 4, 2).
- Install a fluorescent illuminator (fig. 2, 17) on the microscope stand flange (fig. 1, 12). When installing the illuminator, first press the conical surface of the slip on flange to two supports arranged on the right in the stand seat, then clamp the flange with a screw (fig. 1, 28).
- Install the knob (fig. 2, 1) in the position of beam of rays interception with a shutter, having extended it from the casing.
- Place the mercury-filled lamp flashlight (fig. 1, 8) on a table, unscrew the cap clamping screw (fig. 1, 9) and remove the cap (fig. 5, 1).
- Take the mercury-filled lamp from the microscope package, install it into bushings (fig. 5, 2 and 5) on the cap (fig. 5, 1) and fix with screws (fig. 5, 3 and 6).

WARNING! DO NOT TOUCH THE MERCURY-FILLED LAMP BULB! AFTER LAMP INSTALLATION, DEGREASE THE BULB SURFACE WITH ALCOHOL SOLUTION.

- Install the cap (fig. 5, 1) into the mercury-filled lamp flashlight (fig. 1, 8) and fix with a screw (fig. 1, 9).
- Install the mercury-filled lamp flashlight (fig. 1, 8) on the fluorescent illuminator (fig. 2, 15), using the fixator and the adjusting ring (fig. 5, 7 and 8), fix with a bayonet ring located on the illuminator.
- Connect the cable from the flashlight to the slot on the rear surface of the mercury-filled lamp power supply unit (fig. 6).
- Connect the power cord to the power socket on the rear surface of the mercury-filled lamp power supply unit. Make sure that the switch is in "O" position (fig. 1, 18).
- Install the trinocular head (fig. 1, 1) on the fluorescent illuminator flange (fig. 2, 15), fix with a lock knob (fig. 1, 3). Install the luminous flux switch (splitter) (fig. 1, 2) into "Observation only" position.
- Install eyepieces (fig. 2, 5) into eyepiece tubes.
- Install the beam splitting units switching ring (fig. 1, 27) to position No.1.
- Put the stage (fig. 2, 10) down by rotation of the coarse focusing mechanism knob (fig. 2, 13) until stop.
- Install objective lenses (fig. 2, 7) into revolving nosepiece seats in the ascending order of their magnifications.
- Turn the lamp filament adjustment knob (fig. 1, 17) in the direction of brightness reduction until stop.
- The switch (fig. 1, 18) should be installed into "O" position.
- Connect the power cord to the power socket on the rear surface of the stand base (fig. 1, 12).
- Install the UV protection screen (fig. 1, 25) and fix it with screws (fig. 1, 26).

Microscope use

Safety precautions

The microscope may be handled by persons with special medical education. The source of hazard in microscope operation is electric current. Microscope design prevents accidental contact with energized current-conducting parts.

WARNING! REPLACE LAMPS IN FLASHLIGHTS WHEN THE MICROSCOPE AND THE MERCURY-FILLED LAMP POWER SUPPLY UNIT ARE DISCONNECTED FROM THE GRID. TO AVOID HAND SKIN BURN FROM THE LAMP BULB, REPLACE THE LAMP 15–20 MINUTES AFTER DISCONNECTION.

When safety fuses are replaced, it is necessary to install new safety fuses with the same ratings as before.

After the end of operation, the microscope and the mercury-filled lamp power supply unit must be disconnected from the grid.

It is not recommended to leave the appliances connected to the grid unattended.

Perform repair and preventive maintenance only after disconnection of the appliances from the grid.

Observation of objects in transmitted light

Halogen lamp activation and illumination setup

Connect the microscope power cord to AC grid.

Activate the halogen lamp, having set the switch (fig. 1, 18) in "I" position.

Adjust lamp brightness by rotation of the filament adjustment knob (fig. 1, 17).

Image quality in the microscope to a large extent depends on illumination, therefore, illumination setup is an important preparatory operation, which must be performed as follows:

- put the object on the stage (fig. 2, 10) of the microscope;
- activate the objective lens with 4x or 10x magnification into the path of rays (it is recommended to start the process of focusing from low or medium magnification objective lenses with sufficiently large fields and operating distances);
- focus the microscope by rotation of knobs (fig. 1, 14 and 15);
- cover the field diaphragm with a ring (fig. 1, 23), and the aperture diaphragm of the condenser – with a knob (fig. 3, 1);

- observing the object image, focus the condenser, moving it along the height with the knob (fig. 3, 5) for sharp image of the iris field diaphragm;
- if the field diaphragm image is displaced, bring the image into the field center with condenser alignment screws (fig. 3, 2);
- open the field diaphragm with a ring (fig. 1, 23) along the eyepiece field diameter – so that edges of the iris diaphragm are slightly beyond the eyepiece field;
- remove the eyepiece from the right eyepiece tube;
- observing the exit pupil image in the right tube, open the aperture diaphragm of the condenser with the knob (fig. 3, 1) to the exit pupil size. Make sure that the lamp filament image fills the eye. If the image is displaced by the lamp position adjustment knobs (fig. 1, 14 and 15), align the filament image. Using the collector adjustment knob (fig. 1, 16), fill the objective exit pupil with light;
- install the eyepiece in the right eyepiece tube.

Normal operation of the illumination system is provided only when slides with thickness of 1–1.2mm are used.

Microscope focusing for binocular observation

Focus the microscope on the object when observing through a binocular tube as follows:

- put the object on the stage (fig. 2, 10) of the microscope;
- put the objective lens with the necessary magnification into the path of rays;
- rotating the coarse focusing knob (fig. 2, 13), carefully lift the stage to the distance of 0.5mm to the objective lens;
- observing with your right eye into the eyepiece installed into the right eyepiece tube, slowly put the stage down, rotating the coarse focusing knob (fig. 2, 13). As the object contours appear, focus the microscope using a fine focusing knob (fig. 1, 19 or fig. 2, 12);
- observing with your left eye (the right eye is closed) into the eyepiece installed in the left eyepiece tube, get the sharp image of the object by rotation of the diopter mechanism ring (fig. 2, 4). Do not touch the focusing mechanism knobs when doing this;
- set the distance between axes of eyepiece tubes of the binocular head in accordance with the observer's eye base by rotation of the casings with eyepiece tubes in respect to the hinged joint axis so that images of the object in each eyepiece of the head when observing with two eyes are perceived by the observer as one;
- start investigating the slide.

To achieve the best image quality, it is recommended to close the aperture diaphragm of the condenser by 1/3 of the objective exit pupil for each objective lens.

Selection of objective lenses

It is recommended to investigate the object from the objective lens of the lowest magnification, which is used as a search objective lens when selecting a site for more detailed investigation.

After the site for investigation was selected, put its image into the microscope field center. If this operation is performed not carefully enough, the object site of interest for the observer may not get into the field of the stronger objective lens as magnifications are changed.

Then you can go to work with stronger objective lenses, including an objective lens for oil immersion.

Immersion objective lens handling

Handle the immersion objective lens in the premises with temperature from +15 to +25°C (59 to 77°F). Use immersion oil with refraction index $n_D = 1.515$.

WARNING! DO NOT USE SUBSTITUTES INSTEAD OF IMMERSION OIL AS IT MAY SIGNIFICANTLY WORSEN IMAGE QUALITY.

Prior to immersion objective lens handling, set up the microscope as specified in subsections "Halogen lamp activation and illumination setup" and "Selection of objective lenses" and clearly identify the object site for more detailed investigation.

To handle the immersion objective lens:

- lower the stage with a knob (fig. 2, 13);
- apply immersion oil on the object;
- carefully lift the stage using coarse focusing knobs (fig. 2, 13) until the objective lens contacts with the immersion drop on the object;
- observing into the eyepieces and using the fine focusing knob (fig. 1, 19 or fig. 2, 12), get the sharp image of the examined object.

If images of air bubbles that may be contained in the immersion oil layer appear in the eyepiece field in process of focusing, use coarse focusing knobs (fig. 2, 13), put the stage down and repeat focusing.

To investigate the objects, it is necessary to make sure that the iris diaphragm of the objective lens 100x/1.3 is open.

To increase the image contrast, first adjust the condenser aperture diaphragm with a knob (fig. 3, 1), then perform finer adjustment of the contrast with the objective lens iris diaphragm.

Upon completion of operation, remove immersion oil from the frontal objective lens with blotting paper and wipe the contaminated surfaces with cotton wrapped on a stick and slightly soaked in ether or alcohol mix.

Do not push on the frontal objective lens when cleaning.

If the image contrast has lessened or sharpness disappeared as a result of wrong handling of the immersion objective lens, the following is recommended:

- unscrew the objective lens, clean the frontal lens as shown above;
- using oblique light from a table-lamp and a magnifying glass, make sure that there is no dirt, traces of immersion oil, cracks or dents on the frontal lens surface;

- check the microscope illumination setup (the condenser aperture diaphragm shall be open by size of the lens eye or by 2/3 of the eye size).

Observation of objects in fluorescence light

Mercury-filled lamp activation and illumination setup

Connect the mercury-filled lamp power supply unit to the grid. Activate the mercury-filled lamp by putting the power switch in "I" position.

It takes at least 10 minutes for the mercury-filled lamp to reach operating parameters. Normal mode of lamp operation means that amperemeter and voltmeter arrows are in the middle of the scale.

WARNING! DO NOT DEACTIVATE THE MERCURY-FILLED LAMP EARLIER THAN 15 MINUTES AFTER IGNITION! YOU CAN REACTIVATE THE LAMP ONLY 15–20 MINUTES AFTER ITS DEACTIVATION!

Draw "+" sign on a sheet of white paper sized as the stage and put the sheet on the stage. Put the 4x magnification objective lens into the path of rays. Slide the knob (fig. 2, 1) into the casing and put the filter in the path of rays. Using knobs (fig. 1, 5 and 6), open the field and aperture diaphragm. Install the ring (fig. 1, 27) to switch spectral beam splitting units in position No.3 ("B").

Observing into the eyepiece and moving the stage with the coarse focusing knob (fig. 2, 13), get the image of the paper sheet surface. Moving the paper sheet on the stage surface, bring the "+" image into the eyepiece field center. Put the revolver seat free of the objective lens into the path of rays.

Observing from the side (not into the eyepiece) on the paper sheet surface, move the collector with the knob (fig. 1, 7), to get the sharpest image of the mercury-filled lamp discharge arc and its electrodes. Using knobs 10 and 11 (fig. 1), which regulate the mercury-filled lamp position, put the discharge arc image onto "+" sign on the paper sheet surface (in the eyepiece field center). Put the 4x and then 10x magnification objective lens into the path of rays. Observing into the eyepiece and moving the collector with the knob (fig. 1, 7), achieve the most even illumination of the field.

Observation of objects

For investigation in the fluorescence light, the objects are exposed to treatment with special dyes (fluorochromes) having specific spectral characteristics of absorption (excitation) and glow. In accordance with the fluorochrome used to treat the slide, it is necessary to install one of five beam splitting units into the path of rays of the fluorescent illuminator given in table 1 of subsection "Impinging light illumination system".

For example, for rather common treatment of slides using FITC, spectral unit No.3 ("B") is required, for Auramine – unit No.4 ("V"), for stains glowing in the red spectral range, – unit No.2 ("G"). For DAPI and Hoechst stains, unit No.6 ("U") is used; when working with this unit, the filter must be removed from the path of rays with a knob (fig. 2, 1).

Further do the following:

- install the object on the stage (fig. 2, 10) of the microscope;
- activate the objective lens with 10x magnification into the path of rays (it is recommended to start the process of focusing from low or medium magnification objective lenses with sufficiently large fields and operating distances);
- focus the microscope by rotation of knobs (fig. 2, 12 and 13) for the sharp image of the object;
- cover the field and aperture diaphragm with knobs (fig. 1, 5 and 6);
- observing the object image, make sure that the image of the iris field diaphragm is located concentrically to the eyepiece field (if the diaphragm is displaced, align it);
- if the field diaphragm image is displaced, put the image in the center of the field using knobs (fig. 2, 2);
- open the field diaphragm with a knob (fig. 1, 5) along the eyepiece field diameter so that edges of the iris diaphragm are slightly beyond the eyepiece field;
- open the aperture diaphragm with a knob (fig. 1, 6), observing the eyepiece field, make sure that the illumination is quite even, if necessary, adjust the collector focusing with a knob (fig. 1, 7);
- perform focusing on the object for observation with binocular tubes in the same manner as shown in subsection "Microscope focusing for binocular observation" when working in transmitted light;
- start investigating the objects with short breaks in your work. To prevent slide fading, it is necessary to intercept the luminous flux from the lamp with a knob (fig. 2, 1).

Microscope magnification and field diameter on the object

The total microscope magnification Γ in process of visual observation with a binocular head is determined using the following formula:

$$\Gamma = B_{ob} \cdot B_h \cdot \Gamma_{eye}$$

where B_{ob} – linear magnification of the microscope objective lens;

B_h – linear magnification of the head equal to 1.0;

Γ_{eye} – visible magnification of the eyepiece.

The field diameter observed on the object, D_{ob} mm, is determined using the following formula:

$$D_{ob} = \frac{D_{eye}}{B_{ob} \cdot B_h}$$

where D_{eye} – diameter of the eyepiece field limited by the field diaphragm of the eyepiece, in mm.

Possible faults of the microscope and methods of their elimination

External manifestation of the fault	Probable cause	Method of elimination
Cut or uneven illumination	The revolving nosepiece is not in fixation position (the objective lens is not on the microscope axis)	Tighten the revolving nosepiece and put the objective lens into the fixed position, i.e. on the optical axis
	Any objective or eyepiece lens etc. is contaminated	Visually inspect lenses and clean them
	The condenser is not in operating position – too low or warped	Put the condenser in the operating position
There is dust, dirt in the field	Any lens or the stage is contaminated	Remove dirt
Poor quality of object image (low resolution, poor contrast)	There is no cover slip on the object, or its thickness does not comply with the standard	Use the object with a cover slip of standard thickness of 0.17mm
	The object is put with its cover slip down	Turn the object over
	Immersion oil got onto the frontal objective lens. There is dry immersion oil on frontal objective lens 100x \approx /0.17	Clean the immersion oil from the surfaces of frontal objective lenses
	Immersion oil was not applied on the 100x frontal objective lens	Apply oil
	There are bubbles in the immersion oil	Clean the immersion oil from the objective lens, the object, the stage and reapply it
	The aperture diaphragm of the condenser is open or closed too much	Set the necessary diaphragm size
Object images when observed with two eyes in two eyepieces do not match	Eyepiece tubes of the binocular head are not set in accordance with the observer's eye base	Install the binocular head in accordance with the instructions of subsection "Microscope focusing for binocular observation"
When switching the low magnification objective lens to higher magnification objective lens, the objective lens hits the object	The stage with the object is turned over	Install the object with the stage up
	The cover slip is too thick	Use the cover slip of standard thickness
The halogen lamp is not on after activation	The lamp is out. The fuse blew out (safety fuse).	Replace the lamp in accordance with the instructions of subsection "Transmitted light illumination system" Disconnect the microscope from the grid and replace the fuses
The mercury-filled lamp does not go on or went off	The power supply unit is off	Check POWER indication on the power supply unit casing; if it is missing, disconnect from the grid and replace safety fuses with the new ones from the package
	The mercury-filled lamp is installed wrongly	Disconnect the power supply unit from the grid, disconnect the flashlight cord from the unit. Remove the flashlight (after it cooled down), check the proper installation of the lamp according to the instructions of subsection "Impinging light illumination system"
The object fluorescence brightness has reduced greatly	The lamp failed – the bulb has become blurred	Replace the lamp following the instructions of subsection "Impinging light illumination system"

Specifications

Microscope type	biological/optical
Research method	fluorescence, bright field
Magnification, x	40–400
Optics material	optical glass
Eyepiece head	trinocular, inclined at 30°, with luminous flux switching (splitting)
Interpupillary distance, mm	50–75
Eyepiece tube diameter, mm	30
Third vertical eyepiece tube, mm	23.2
Eyepieces	wide field WF 10x/22mm with eyecups (2pcs.)
Objective lenses	infinite semi-apochromatic luminescent (fluorescent), objective lenses: 4x, 10x, 20x, 40x
Revolving nosepiece	6 objectives
Stage	mechanical double-layer, 180x160mm, with a mechanical scale
Stage moving range, mm	85x50
Condenser	detachable Abbe condenser N.A. 1.25 with an iris diaphragm and a filter holder
Diaphragm	iris, field
Focusing	coaxial, coarse and fine focusing; fine focusing scale: 0.002mm
Body material	metal
Illumination	lower (12V/30W halogen lamp) with brightness adjustment
Fluorescent module	filters "G", "B", "BV", "V", "U"; mercury-filled lamp (100W) with external power supply unit; UV protection screen
Power supply	AC adapter 100–220V/50–60Hz

Levenhuk reserves the right to make changes to the product range and specifications without prior notice.

CAUTION! PLEASE REMEMBER THAT MAINS VOLTAGE IN MOST EUROPEAN COUNTRIES IS 220–240V. IF YOU WANT TO USE YOUR DEVICE IN A COUNTRY WITH A DIFFERENT MAINS VOLTAGE STANDARD, REMEMBER THAT USE OF A CONVERTER IS ABSOLUTELY NECESSARY. THE MICROSCOPE MUST BE GROUNDED. MAKE SURE THAT THE MAIN VOLTAGE MATCHES THE VOLTAGE INDICATED ON THE MICROSCOPE BODY.

Care and maintenance

- **Never, under any circumstances, look directly at the Sun, another bright source of light or at a laser through this device, as this may cause PERMANENT RETINAL DAMAGE and may lead to BLINDNESS.**
- Take necessary precautions when using the device with children or others who have not read or who do not fully understand these instructions.
- After unpacking your microscope and before using it for the first time check for integrity and durability of every component and connection.
- Do not try to disassemble the device on your own for any reason. For repairs and cleaning of any kind, please contact your local specialized service center.
- Protect the device from sudden impact and excessive mechanical force. Do not apply excessive pressure when adjusting focus. Do not overtighten the locking screws.
- Do not touch the optical surfaces with your fingers. To clean the device exterior, use only special cleaning wipes and special optics cleaning tools from Levenhuk. Do not use any corrosive or acetone-based fluids to clean the optics.
- Abrasive particles, such as sand, should not be wiped off lenses, but instead blown off or brushed away with a soft brush.
- Do not use the device for lengthy periods of time, or leave it unattended in direct sunlight. Keep the device away from water and high humidity.
- Be careful during your observations, always replace the dust cover after you are finished with observations to protect the device from dust and stains.
- If you are not using your microscope for extended periods of time, store the objective lenses and eyepieces separately from the microscope.
- Store the device in a dry, cool place away from hazardous acids and other chemicals, away from heaters, open fire and other sources of high temperatures.
- When using the microscope, try not to use it near flammable materials or substances (benzene, paper, cardboard, plastic, etc.), as the base may heat up during use, and might become a fire hazard.
- Always unplug the microscope from a power source before opening the base or changing the illumination lamp. Regardless of the lamp type (halogen or incandescent), give it some time to cool down before trying to change it, and always change it to a lamp of the same type.
- Always use the power supply with the proper voltage, i.e. indicated in the specifications of your new microscope. Plugging the instrument into a different power outlet may damage the electric circuitry of the microscope, burn out the lamp, or even cause a short circuit.
- **Seek medical advice immediately if a small part or a battery is swallowed.**

Levenhuk International Lifetime Warranty

All Levenhuk telescopes, microscopes, binoculars, and other optical products, except for their accessories, carry a **lifetime warranty** against defects in materials and workmanship. A lifetime warranty is a guarantee on the lifetime of the product on the market. All Levenhuk accessories are warranted to be free of defects in materials and workmanship for **six months** from the purchase date. The warranty entitles you to the free repair or replacement of the Levenhuk product in any country where a Levenhuk office is located if all the warranty conditions are met.

For further details, please visit: www.levenhuk.com/warranty

If warranty problems arise, or if you need assistance in using your product, contact the local Levenhuk branch.

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
The original Levenhuk cleaning
accessories

Levenhuk Cleaning Pen LP10

Removes dust with a brush
The soft tip is treated with a special cleaning fluid that removes greasy stains
Does not damage optical coatings of the lenses
Leaves no smudges or stains

levenhuk.com

Levenhuk Inc. (USA): 928 E 124th Ave. Ste D, Tampa, FL 33612,
USA, +1 813 468-3001, contact_us@levenhuk.com
Levenhuk Optics s.r.o. (Europe): V Chotejně 700/7, 102 00 Prague 102,
Czech Republic, +420 737-004-919, sales-info@levenhuk.cz
Levenhuk® is a registered trademark of Levenhuk, Inc.
© 2006–2023 Levenhuk, Inc. All rights reserved.
20230529

levenhuk
Zoom&Joy